


10197

# Rubber Clip x 4 for 1017 and 1018, 120 mm, Orange



4 hooks and 8 securing pins for wall brackets 1017x and 1018x.  
Rubber clips can hold products with diameters of 25–35 mm.

# Technical Data

<b>Item Number</b>	10197
<b>Material</b>	TPE Rubber
<b>Complies with (EC) 1935/2004 on food contact materials<sup>1</sup></b>	No
<b>Produced according to EU Regulation 2023/2006/EC of Good Manufacturing Practice</b>	Yes
<b>FDA compliant raw material (CFR 21)</b>	Yes
<b>Use of phthalates and bisphenol A</b>	No
<b>Box Quantity</b>	5 Pcs.
<b>Quantity per Pallet (80 x 120 x 200 cm)</b>	3200 Pcs.
<b>Quantity Per Layer (Pallet)</b>	0 Pcs.
<b>Box Length</b>	300 mm
<b>Box Width</b>	80 mm
<b>Box Height</b>	20 mm
<b>Length/Depth</b>	120 mm
<b>Width</b>	80 mm
<b>Height</b>	20 mm
<b>Gross Weight</b>	0,06 kg
<b>Net Weight</b>	0,05 kg
<b>Cubik metre</b>	0,0002 M3
<b>Recommended sterilisation temperature (Autoclave)</b>	121 °C
<b>Max. cleaning temperature (Dishwasher)</b>	93 °C
<b>Max usage temperature (non food contact)</b>	80 °C
<b>Min. usage temperature<sup>3</sup></b>	-20 °C
<b>Min. pH-value in usage concentration</b>	2 pH
<b>Max. pH-value in Usage Concentration</b>	10,5 pH
<b>GTIN-13 Number</b>	5705020101977
<b>GTIN-14 Number (Box quantity)</b>	15705020101974
<b>Customs Tariff No.</b>	39241000
<b>Country of origin</b>	Denmark

New equipment should be cleaned, disinfected, sterilised and any labels removed, as appropriate to its intended use, e.g. high risk vs. low risk food production areas, general hospital areas vs. intensive care units, before use.